 CEDARS-SINAI BIOMANUFACTURING CENTER	CEDARS SINAI BIOMANUFACTURING CENTER INDUCED PLURIPOTENT STEM CELL CORE THE DAVID and JANET POLAK FOUNDATION STEM CELL CORE LABORATORY	iPSC PASSAGING PROTOCOL WITH PICK-TO-KEEP	
		SOP NUMBER: SOP-iPSC-010	Version: A

1. PURPOSE

To describe the procedure for passaging of iPSCs maintained on Matrigel using Pick-to-Keep (PTK)

2. SUPPLIES

Complete mTeSR Plus Medium (Basal medium + 5x Supplement) (StemCell Technologies, Cat # [05825](#))

Matrigel Coated TC dish (Prepared as described in SOP-iPSC-002)

BD Lo-Dose™ U-100 Insulin Syringes (FisherSci, Cat [#329461](#))

5ml and 10ml sterile serological pipettes

Filtered 1000uL Pipette tips

Unfiltered 200uL Pipette tips

Biohazardous Sharps Waste Container (inside BSC)


3. PROCEDURE

NOTE: You must have a prepared Matrigel coated plate before starting this protocol. If you are using a Matrigel coated plate that has been stored at 4°C, **the plate must be allowed to equilibrate to room temperature for 1 hour prior to starting.**

- 3.1. Aspirate Matrigel from the prepared Matrigel coated tissue culture plate and add 1 ml of complete mTeSR media to each well
- 3.2. You should have a total volume of at least 2.5mls per well after cells have been added
- 3.3. Aspirate spent media.
- 3.4. Add 1ml of mTeSR per wells needed for PTK.
- 3.5. Working under microscope inside the BSC at 4X objective, use an insulin syringe needle to cut iPSC colonies into small clumps (see Figure A)

NOTE: Uncap the syringe one time and discard the cap right away. Avoid injury by never re-capping a needle

- 3.5.1. Take great care to use the lightest pressure possible in order to avoid scratching through the coating substrate and plasticware, as the insulin syringe is incredibly sharp
- 3.5.2. It is recommended to select at least 10 colonies per well of a 6-well plate to obtain sufficient adhesion and proliferation the following week

 CEDARS-SINAI BIOMANUFACTURING CENTER	CEDARS SINAI BIOMANUFACTURING CENTER INDUCED PLURIPOTENT STEM CELL CORE THE DAVID and JANET POLAK FOUNDATION STEM CELL CORE LABORATORY		iPSC PASSAGING PROTOCOL WITH PICK-TO-KEEP
	SOP NUMBER: SOP-iPSC-010	Version: A	

- 3.6. Carefully dispose the used syringe into a biohazardous sharp waste inside the BSC
- 3.7. Place a P1000 micropipette tip inside a P200 tip.
- 3.8. Use the P1000+P200 tip to gently nudge or lift the cut colony pieces off the plate similar to a fine cleaning motion (see Figure B)
- 3.9. Add 1 ml of mTesR to the well with picked clumps and gently rinse the well to collect all the pieces
- 3.10. Tilt the PTK plate and gently resuspend to rinse the well and pool the clumps at the bottom of the well
- 3.11. Transfer floating clumps to new plate.
- 3.12. Place the plate in the 37°C incubator with 5% CO₂ and gently rock the plate back and forth and side-to-side to ensure even distribution of the colonies throughout the well.
- 3.13. Do not move the plate for 24 hours
- 3.14. The next day, check for colony attachment and proliferation

Figure A



Figure B

